

MEDIC-80

DISK INSTRUCTIONS:

This program is supplied on a FORMATTED disk. It will run best under TRSDOS 1.3.

To Start From TRSDOS Ready*.....
Type DO MEDIC

To Start From BASIC Ready
>
Type RUN"MEDONE"

TAPE INSTRUCTIONS:

If you order the Model III version it is recorded at 1500 baud and will load under the CLOAD command.

If you order the Model I version it is recorded at 500 baud and also will load under the CLOAD command.

After loading type RUN and follow the instructions on the screen. The program is in two sections recorded on each side of the cassette.

PROGRAM COMMENTS

You are ultimately responsible for your medical health. This program is limited in scope to help deal with 90% of the medical questions that a person will ask. If you suspect a major medical problem please check with your doctor. This program suggests solutions only to problems that can usually be dealt with in the home. When in doubt about your responses it will suggest a consult with a doctor. No other warranties are express or implied. And no liability will be accepted by the selling agent for ability or inability to use this program.

It is completely MENU driven (just select either number or Y or N) as the program requests. Emergencies and skin problems are covered in the first module with common symptoms being covered in the second module. The program will pause while the second or first module loads into the computer. People using the tape version please follow the on screen directions!

The program can also be used to test your knowledge in any areas. For example explore as you like the Emergency section, you will most likely find some new and interesting information located there! Remember that at any level in the program you can return to the TOP MENU for that module by press "Q" for quit and return top.

* TRSDOS Copyright of Tandy Inc.